

BEYOND™ WhiteSpa SELECT

Cosmetic Teeth Whitening System

beyond™ Dental & Health

BEYOND™ USA HEADQUARTERS
12902 Devonage Dr. Suite 500
Dallas, TX 75243 USA
Tel: +1-214-271-0100
Fax: +1-214-271-4288
E-mail: usa@beyonddental.com

BEYOND™ EUROPEA HEADQUARTERS
Jugend 17-494-20
12107 Berlin, Germany
Tel: +49-30-3060-800
Fax: +49-30-3060-801
E-mail: europa@beyonddental.com

BEYOND™ ASIA-PACIFIC HEADQUARTERS
No. 85, Chin Avenue, #12-02/0211
Singapore
Tel: +65-6336-1000
Fax: +65-6336-1700
E-mail: asia@beyonddental.com
<http://www.beyonddental.com>

©2008 Beyond Dental & Health. All rights reserved. • 1034 p.m. 03/08

Q: Who can undergo the BEYOND™ WhiteSpa Select treatment?

A: BEYOND™ WhiteSpa Select produces dramatic whitening results for anyone who wants to whiten his or her teeth. This includes people with teeth that have been stained by smoking or by substances such as coffee or tea. It also whitens teeth stained by tetracycline, specified by fluoride, or hereditary discoloration. Pregnant or lactating women or clients under 16 years of age should not receive the BEYOND™ WhiteSpa Select treatment.

Q: Is it safe?

A: Yes, the BEYOND™ WhiteSpa Select treatment is safe. We use a gentle, over-the-counter strength peroxide-based whitening gel the ONLY proven safe whitening agent combined with a tailored light source that removes all harmful UV light and heat before delivering a concentrated beam of light to accelerate the whitening process.

Q: Why is the BEYOND™ WhiteSpa Select treatment one hour?

A: The BEYOND™ WhiteSpa Select treatment is designed to produce long-lasting, dramatic whitening results. Effective teeth-whitening requires the combination of a powerful, tailored light source with a specially formulated, quick-reacting, hydrogen peroxide whitening gel exposed to the teeth for a minimum of 60 minutes. Cheaper systems that combine weak LED lights with a slow-reacting, carbamide peroxide whitening gel and advertise shorter treatment times will NOT provide the same level of quality, long-lasting whitening as the BEYOND™ WhiteSpa Select.

Q: Will it whiten my existing crowns or veneers?

A: Whitening treatments can only whiten natural teeth and cannot change the color of porcelain crowns, veneers, or fillings.

Q: How white will my teeth be after the treatment?

A: Whitening results can vary by individual depending on the severity and nature of your stains. Your BEYOND™ WhiteSpa Select treatment will be tailored to produce the best whitening results possible. Talking to your aesthetician about your expectations from the treatment will assist them in designing an appropriate whitening schedule to meet your individual needs.

Q: How long will the whitening results last?

A: This varies from person to person and depends on your diet and lifestyle, but the effects usually last for about two years. Periodic touch-up treatments are recommended for clients who consume coffee or other staining foods and drinks frequently or for those who use tobacco products.

Q: Why is BEYOND™ WhiteSpa Select the best?

A: BEYOND™ WhiteSpa Select is the original premium cosmetic tooth whitening system. Containing a powerful halogen light source (proven to be more effective than LED lights) and an advanced light filtration system, BEYOND™ WhiteSpa Select emits an intense, pure beam of light to effectively accelerate the whitening process. Designed exclusively for the beauty industry, no other cosmetic whitening system produces the same level of whitening results while ensuring your comfort and safety.

Your Brightest
SMILE in just 1 hour!

BEYOND™
WhiteSpa
SELECT
Cosmetic Teeth Whitening System

BEYOND™ WhiteSpa SELECT

Cosmetic Teeth Whitening System

No more messy whitening trays or strips and uncomfortable dental-office treatments!

Enjoy a premium whitening experience from the comfort of your favorite spa with BEYOND™ WhiteSpa Select and discover a reason to smile with confidence again!

» Ask your aesthetician about the BEYOND™ WhiteSpa professional tooth whitening treatment today!

Introducing your solution for a bright, white, confident smile:

BEYOND™ WhiteSpa Select is in the business of making people smile with confidence. Whether you have a big job interview, graduation, wedding, or event coming up – or if you just want to put your best face forward everyday – having shiny, white teeth will restore your confidence in letting the world see your beautiful smile.

BEYOND™ WhiteSpa Select Professional Cosmetic Tooth Whitening

Select the only true power whitening system available through the beauty salon or spa – and select the best!



before

after

BEYOND™ WhiteSpa Smile with Confidence!

As a global leader in the design and manufacture of professional teeth whitening products, BEYOND™ is dedicated to developing the safest and most advanced technologies for cosmetic teeth whitening.

The BEYOND™ WhiteSpa Select system produces a concentrated beam of high-intensity blue light to accelerate the whitening process. An advanced filtration system removes harmful ultraviolet light and heat, while BEYOND™ WhiteSpa Select's specially formulated whitening gel works to whiten teeth, restoring them to a shiny white.

In just three, 20 minute cycles, stains that have accumulated over the years are rendered neutral in color and a beautiful, shiny, white smile is revealed.

- ✓ Gentle whitening gel formula reduces risk of sensitivity!
- ✓ Powerful, halogen light source produces optimal whitening results!
- ✓ Zero exposure to harmful UV light or heat!

